



mc2



**MOLINO DE MARTILLOS *MM-100***



### MOLINO DE MARTILLOS *MM-100*

- Martillo y tolva inoxidable y fácilmente desmontables.
- Cribas intercambiables de 4,5 y 5,5 mm.
- Bandeja receptora de pasta.
- Robusto chasis de acero.
- Carenado de poliéster reforzado con fibra de vidrio.
- Motor eléctrico monofásico de 1,5 kW.

Tiempo requerido para tratar una muestra:  
Un minuto.

El Molino de Martillos *MM-100* puede funcionar en combinación con el resto de los elementos del sistema **Abencor**® o de forma individual, bastando disponer para su alimentación de un enchufe normal de corriente alterna de 220-230 V 16 A con toma de tierra incorporada.

### HAMMER MILL *MM-100*

- Easily detachable stainless steel hammer and load funnel.
- Interchangeable screens 4.5 and 5.5 mm.
- Paste reception tray.
- Rugged steel frame.
- Fibre glass reinforced polyester cover.
- Monophase electric 1.5 kW, 220-230 V.

Time required to treat one sample:  
One minute.

The Hammer Mill *MM-100* may work alone or integrated in the **Abencor**® system. A standard electric socket AC 220-230 V 16 A, 50 Hz, earthed, is all you need. Special drives for different voltage and frequency are available under demand.

---

#### MC2 INGENIERÍA Y SISTEMAS, S. L.

Diego Angulo Iñiguez, 14  
Tel. 954 987 388 Fax 954 536 832  
E-mail: emecedos@interbook.net  
41018 SEVILLA -España (Spain).